

Isolated Tissue Studies

PowerLab® and Multiple Chamber Organ Bath Systems



ADInstruments PowerLab systems are widely used to record and analyze data from experiments using isolated tissues. These tissues typically include isolated muscle, arterial rings or strips, uterine muscle, vas deferens, ileum, colon, isolated atrium, ventricle or diaphragm. Our range of multi-chamber organ baths, bridge amplifiers, stimulating electrodes and isometric and isotonic transducers are designed to complement PowerLab systems in this area of research.

Two, four or eight chamber organ bath configurations are available. Each organ bath can include tissue chambers with 5, 10, 25 or 50 mL volume. Pre-heating reservoir coils are located in a compact temperature-controlled water bath constructed from thick Perspex®. This ensures that the temperature within the glass tissue chambers remains constant. Tissue holders, micropositioners, water pump and temperature controller are included.

The solid design, combined with the volume of water, minimizes any signal artifact due to mechanical vibration. Each tissue chamber has an adjustable gas diffuser. Electrovalves enable simultaneous or individual emptying and refilling of tissue chambers at the press of a button.

To complement the hardware, ADInstruments LabChart® software (included with the PowerLab) and analysis modules such as Dose Response, enable laboratories to conduct studies and analysis efficiently to allow faster throughput of experiments.

Features & Benefits

- Complete systems ready for experiments at affordable prices
- Easy filling and emptying of tissue chambers at the push of a button
- Immersed chambers and reservoir coils ensure temperature stability
- Internal heater and thermostat provide stable temperature control
- Compact, robust design with no external glassware
- Quick to connect and easy to use




Data Acquisition & Analysis

PowerLab data acquisition systems (comprising a PowerLab recording unit with LabChart software) are ideal for monitoring and recording signals from isolated tissue experiments. PowerLab systems offer seamless integration with bridge amplifiers, transducers, organ baths and most stimulators. Our bridge amplifiers provide software-controlled zeroing, range and filter selection. You can easily calibrate transducers and recording channels into appropriate units such as grams or Newtons.



The flexible LabChart display includes a split screen feature to view and analyze previously recorded data while still recording new data. The powerful Data Pad facility allows easy extraction of experimental results for analyses, with OLE capabilities built in.

LabChart Software Features

- Control of stimulation frequency and timing of pulse trains
- Online and offline calculations
- Timed Add to Data Pad for fast and easy extraction of experimental results
- Easy calibration (using Units Conversion)
- Comment annotation during or after recording
- Automation of experimental procedures
- Software selectable sampling rates, range and filter settings
- Additional advanced calculations and displays using the LabChart Dose Response Module
- Easy extraction of data into a spreadsheet or other programs, including GraphPad Prism



Above: Response of spontaneously beating guinea pig atria to increasing concentrations of isoproterenol recorded in LabChart.



Above: Relaxation of precontracted rat uterine muscle preparation, with spontaneous contractile activity at left. KCl (40 mM) added at comment #4 precontracts the muscle. Increasing concentrations of salbutamol (10^{-10} M to 10^{-5} M) added at the other comments show a cumulative effect.

The Dose Response Results View displays the measured points and calculated dose response curves. The right hand panel displays the best fit values.

Rat Uterus.adich: Dose Response Table View		
	A	B
	Salbutamol	
		<Run 2 Label>
Dose	Ch1: Tension	<Select Channel>
1	0.1	9.815
2	0.3	9.772
3	1	9.597
4	3	9.106
5	10	9.596
6	30	9.109
7	100	7.842
8	300	4.029
9	1000	3.638
10	3000	3.470
11	10000	3.459

Using the Dose Response Module the Table View displays the dose response data

Product Selection

The following products are recommended for isolated tissue studies. Individual components may be selected or complete organ bath systems can be chosen to reduce cost and provide a ready-to-use solution. See back page for ordering information.



Organ Baths

Compact organ baths come in two, four and eight chamber models. All models feature electrovalves for filling and emptying of tissue chambers individually or simultaneously, at the push of a button. Tissue chambers are available in 5, 10, 25 or 50 mL volumes. Pre-heating reservoir coils, gas diffusers, tissue holders, micropositioners, water pump and thermostat controller are all included.



PowerLab Systems

We recommend eight and sixteen channel PowerLab 35 series models for use with multiple-chamber organ baths. These systems connect via USB to Windows® and Mac® OS computers. The PowerLab data acquisition units feature analog outputs that can be used to control an external stimulator. A trigger input provides remote starting capability. Fast data acquisition, 16-bit resolution, real-time data display and powerful analysis features ensure accuracy and efficiency.



Bridge Amps

ADInstruments Quad and Octal Bridge Amps are optimized for use with PowerLab systems and multi-chamber organ baths. They include fast software or switch initiated zeroing. Each input provides a software-controlled, low-noise bridge amplifier with software-selectable input ranges and filtering. These models are suitable for use with full-bridge strain-gauge transducers.

Transducers

Isometric transducers are available in regular and high sensitivity models. If you want to measure displacement against a constant preload, we also provide isotonic transducers. Micropositioners hold the transducers above the tissue chambers to enable fine adjustment of tissue pretension. ADInstruments supplied transducers plug directly into ADInstruments Bridge Amps for immediate use.



Stimulating Electrodes

We supply single ring, double ring or pole stimulating electrodes in platinum or stainless steel. A choice of internal diameters is available for ring electrodes. You will require a stimulator for use with stimulating electrodes. Your ADInstruments representative can provide information on a suitable stimulator configuration.



Ordering Information

Organ Bath Systems

PL3508B5/C* Organ Bath System	PL3508B6/C* Organ Bath System	PL3516B7/C* Organ Bath System
1 x PL3508/P PowerLab 8/35 includes LabChart Pro software	1 x PL3508/P PowerLab 8/35 includes LabChart Pro software	1 x PL3516/P PowerLab 16/35 includes LabChart Pro software
1 x FE224 Quad Bridge Amp	1 x FE228 Octal Bridge Amp	2 x FE228 Octal Bridge Amps
4 x MLT0420 Force Transducer (20 g)	8 x MLT0420 Force Transducer (20 g)	16 x MLT0420 Force Transducer (20 g)
1 x ML0146/C 4 Chamber Organ Bath	1 x ML0186/C 8 Chamber Organ Bath	2 x ML0186/C 8 Chamber Organ Bath

* Need to specify Chamber size 5, 10, 25 or 50 mL.

Individual Items

Organ Baths**	
ML0126/C Two Chamber Organ Bath and Thermostat Controller	ML0186/C Eight Chamber Organ Bath and Thermostat Controller
ML0146/C Four Chamber Organ Bath and Thermostat Controller	

**Organ baths include: tissue chambers, pre-heating reservoir coils, gas diffusers, tissue holders, micropositioners, water pump and thermostat controller. Tissue chamber size C = 5, 10, 25 or 50 mL, must be selected on ordering.

Transducers			
MLT0420 Force Transducer (20 g)	Operating range: ~0 g up to 20 g	Sensitivity: 87.5 μ V/V/g	
MLT0402 Force Transducer (2 g)	Operating range: ~0 g up to 2 g	Sensitivity: 325 μ V/V/g	
MLT0201 Force Transducer	Operating range: 5 mg to 25 g wt	Sensitivity: 30 mV/V (full range)	
MLT0202 Sensitive Isometric Transducer	Operating range: 0 mg to 25 g wt	Sensitivity: 100 mV/V (full range)	
MLT0015 Isotonic Transducer	Displacement range: $\pm 15^\circ$ (0 to 48 mm)	Sensitivity: 1 mV/V/ $^\circ$ (625 μ V/V/mm)	
MLT7006 Hall Effect Isotonic Transducer	Displacement range: $\pm 15^\circ$ (0 to 60 mm)	Sensitivity: 90 μ V/V/ $^\circ$ (45 μ V/V/mm)	

Stimulating Electrodes Specify internal diameter for ring electrodes: ID = 6.5, 8 or 10 mm		
Type and Description	Platinum	Stainless Steel
Single Ring (6.5, 8 or 10 mm inner diameter) 20 mm between ring and tissue holder	MLA0301/ID	MLA0305/ID
Double Ring (6.5, 8 or 10 mm inner diameter) 20 mm between rings, 6 mm above tissue holder	MLA0302/ID	MLA0306/ID
3 mm Length Pole, 20 mm from the tissue holder	MLA0303	MLA0307
Double Ring (6 mm inner diameter) 15 mm between rings, 12 mm from tissue holder	MLA0304	MLA0308

Software			
MLS060/7 LabChart	MLS330/7 GLP Client and MLS335 GLP Server		
MLS260/7 LabChart Pro	(Includes the modules listed below. Modules are also available for individual purchase.)		
MLS390/7 Dose Response (Win)	MLS310/7 Heart Rate Variability (Win and Mac)	MLS340/7 Cardiac Output (Win)	
MLS065/7 DMT Normalization (Win and Mac)	MLS240/7 Metabolic (Win and Mac)	MLS320/7 Video Capture (Win and Mac)	
MLS370/7 Blood Pressure (Win)	MLS062/7 Spike Histogram (Win and Mac)	MLS395/7 Circadian Analysis (Win)	
MLS360/7 ECG Analysis (Win)	MLS380/7 Peak Analysis (Win)	MLS375/7 PV Loop (Win)	

 Share your data with colleagues. Free LabChart Reader – download to view and analyze LabChart data.

PowerLab, MacLab, LabChart, LabTutor and LabAuthor are registered trademarks and Chart and Scope are trademarks of ADInstruments Pty Ltd. All other trademarks are the property of their respective owners ISTPL09/11

PowerLab systems and signal conditioners meet the European EMC directive. ADInstruments signal conditioners for human use are approved to the IEC60601-1 patient safety standard and meet the CSA C22.2 No. 601.1-M90 and UL Std No. 2601-1 safety of medical electrical equipment standards.



ADINSTRUMENTS.com

ISO 9001:2008 Certified Quality Management System 

North America

Tel: +1 888 965 6040
Fax: +1 719 576 3971
info.na@adinstruments.com

United Kingdom

Tel: +44 1865 332 050
Fax: +44 1865 332 051
info.uk@adinstruments.com

Germany

Tel: +49 6226 970105
Fax: +49 6226 970106
info.de@adinstruments.com

North Asia

Tel: +86 21 5830 5639
Fax: +86 21 5830 5640
info.cn@adinstruments.com

South East Asia

Tel: +60 3 8024 5296
Fax: +60 3 8023 6307
info.sea@adinstruments.com

Japan

Tel: +81 52 932 6462
Fax: +81 52 932 6755
info.jp@adinstruments.com

South America

Tel: +56 2 356 6749
Fax: +56 2 356 6786
info.cl@adinstruments.com

Brazil

Tel: +55 11 3266 2393
Fax: +55 11 3266 2392
info.br@adinstruments.com

South Asia

IN Tel: +91 11 4306 5615
PK Tel: +92 21 3489 2518
info.in@adinstruments.com

Australia

Tel: +61 2 8818 3400
Fax: +61 2 8818 3499
info.au@adinstruments.com

New Zealand

Tel: +64 3 477 4646
Fax: +64 3 477 4346
info.nz@adinstruments.com

International

Tel: +61 2 8818 3400
Fax: +61 2 8818 3499
info.au@adinstruments.com